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
2024 Training-Transfer Program for Medical Technicians/Parasite Specialists

Diagnostic Training Transfer Program for Common & Opportunistic Intestinal Parasites Using An Integrated Staining Method

Safranin O Staining


Chromotrope 2R staining

1. Use Chang's stool collection container to collect peanut-sized amount of stool.
2. Add 5ml 10% formalin solution, mix well, filter the solution and collect the filtrate in 15ml tube.
3. Add 5ml 10% formalin. Mix well and centrifuge at 2000rpm for 5 minutes.
4. Take one drop of the precipitate and make an air-dried smear slide.

5. Fix by methanol for 5 minutes.


6. Stain the slide with boiled Safranin O solution for 1 minute and wash by distilled water.



7. Stain the background with methylene blue solution for 2 minutes, wash by water and air-dry the slide. Mount with a coverslip using desired mounting media.

8. Examine under the 100X oil immersion objective lens. Cyclospora oocysts and Cryptosporidium parvum will be stained orange stably under Safranin O staining.

5. Fix by methanol for 1 minute.


6. Place in chromotrope stain for 90 minutes.

7. Destain in acid alcohol for only 1 to 3 seconds.


8. Rinse in 95% ethanol by dipping.

9. Place two changes of 100% ethanol for 3 minutes each.

10. Place two changes of xylene or xylene substitute for 10 minutes each.


11. Drain slide and mount with coverslip using mounting media (e.g., permount). Examine smear after drying using at least 100x objective oil immersion or higher. Examine at least 200 to 300 oil immersion fields.

Microscopic examination of parasite ova in feces by Merthiolate-Iodine-Formaldehyde (MIF) method

Materials

- A. Chang's stool container
B. PP tube & Dropper
C. MIF Reagent
D. Slides

Procedures

1. Use stir bar to collect stool specimen (thumb size)
2. Add 5 to 7 ml of MIF solution
(MF stock : Lugol solution = 15 : 1)
3. Mix well and elute the solution into PP tube
4. Stain the feces overnight
(at least 2 hours is recommended)
5. Examine the parasites ova by microscopy.
For helminths: 100X is suggested for screen; 400X is suggested for confirmation
For intestinal protozoa: 400X is suggested for screen and confirmation

Positive Results

- Please refer to “Color Atlas of Common Parasites in the Kingdom of Eswatini”.

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Questionnaire for Satisfaction Investigation (2024)

1. Do you think this [Common & Opportunistic Intestinal Parasites by An Integrated Staining training-transfer program] is useful to your future lab work?
☐Excellent (10) ☐Good(8) ☐Fair(6) ☐No(4)
2. How do you think about the Topics in this workshop?
☐Excellent (10) ☐Good(8) ☐Fair(6) ☐No(4)
3. Do you think the protocol of the Integrated Staining method is useful to make you be familiar with these new diagnostic systems?
☐Excellent (10) ☐Good(8) ☐Fair(6) ☐No(4)
4. Do you think the [Common & Opportunistic Intestinal Parasites by An Integrated Staining Identification skill] (training course) is useful to make you familiar with the identification skills of i.e. Ascaris, Hookworm, Entamoeba, Cryptosporidium and microsporidia?
☐Excellent (10) ☐Good(8) ☐Fair(6) ☐No(4)
5. Do you think the lab practice is useful to make you more familiar with the performing procedure and identification skills on common & opportunistic intestinal parasites?
☐Excellent (10) ☐Good(8) ☐Fair(6) ☐No(4)
6. Do you think the booklet for this workshop is very nice to be understandable and readable?
☐Excellent (10) ☐Good(8) ☐Fair(6) ☐No(4)
7. Will you consider inviting Prof. Fan with his team to continue providing the [Diagnostic-Training Transfer Program] workshop next year?
☐Absolutely (10) ☐Sure(8) ☐Maybe(6) ☐No(4)
8. Overall, what's your evaluation of this 2-days/1 day workshop?
☐Excellent (10) ☐Good(8) ☐Fair(6) ☐No(4)


Thank you very much for your participation and time!

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
2024 Diagnostic Training Transfer Program for Common & Opportunistic Intestinal Parasites Using An Integrated Staining Method Satisfaction Survey


1.Do you think this [Common & Opportunistic Intestinal Parasites by An Integrated Staining training-transfer program] is useful to your future lab work?


2.How do you think about the Topics in this workshop?


3.Do you think the protocol of the Integrated Staining method is useful to make you be familiar with these new diagnostic systems?


4.Do you think the [Common & Opportunistic Intestinal Parasites by An Integrated Staining Identification skill] (training course) is useful to make you familiar with the identification skills of i.e. Ascaris, Hookworm, Entamoeba, Cryptosporidium and micro


5.Do you think the lab practice is useful to make you more familiar with the performing procedure and identification skills on common & opportunistic intestinal parasites?


6.Do you think the booklet for this workshop is very nice to be understandable and readable?

7.Will you consider inviting Prof. Fan with his team to continue providing the [Diagnostic-Training Transfer Program] workshop next year?

8.Overall, what's your evaluation of this 2-days/1 day workshop?

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Mrs. Sindisiwe Dlamini, Director of Swaziland Health Laboratory Services

Mrs. Gugu Maphalala, Technical Director of Swaziland National Blood Transfusion Service

Mrs. Quinton Dlamini, Director of National Malaria & NTDs Program

Mrs. Fortunate Lushaba, Technician of Swaziland Health Laboratory Services

Mr. Bongani Dlamini, Associate of Mbabane Government Hospital

Mr. Phumlani Motsa, Technician of National Bilharzia (Schistosomiasis) Worm Control Program

Mbabane Government Hospital, Hhohho, Swaziland

Raleigh Fitkin Memorial Hospital, Manzini, Swaziland

Hlatikhulu Government Hospital, Shiselweni, Swaziland

Piggs Peak Government Hospital, Hhohho, Swaziland

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1. Morton S, Pencheon D, Squires N. Sustainable Development Goals (SDGs), and their implementation: A national global framework for health, development and equity needs a systems approach at every level. *Br Med Bull.* 124(1):81-90, 2017.
2. Quintero-Betancourt W, Peele ER, Rose JB. *Cryptosporidium parvum* and *Cyclospora cayetanensis*: a review of laboratory methods for detection of these waterborne parasites. *J Microbiol Methods.* 49(3):209-24, 2002.
3. Field AS, Milner DA Jr. [Intestinal microsporidiosis](#). *Clin Lab Med.* 35(2):445-59, 2015.
4. Araniy F, Faisal, Abdullah A, Bokhari. *Microsporidium*. In: StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2020 Jan. 2020 Aug 13.
5. Chukwuma C Sr. *Microsporidium* in AIDS patients: a perspective. *East Afr Med J.* 73(1):72-5, 1996.
6. Weber R, Bryan RT. Microsporidial infections in immunodeficient and immunocompetent patients. *Clin Infect Dis.* 19(3):517-21, 1994.
7. Collins PA, Wright MS. Emerging intestinal protozoa: a diagnostic dilemma. *Clin Lab Sci.* 10(5):273-8, 1997.
8. Stentiford GD, Ross SH, Kerr R, Bass D, Bateman KS. *Paradoxium irvingi* n.gen. n.sp. (Microsporidia) infecting the musculature of European pink shrimp *Pandalus montagui*. *J Invertebr Pathol.* 130:1-8, 2015.
9. Galvan AL, Magnet A, Izquierdo F, Fenoy S, Rueda C, Vadillo CF, Henriques-Gil N, del Aguila C. Molecular Characterization of Human-Pathogenic Microsporidia and *Cyclospora cayetanensis* Isolated from Various Water Sources in Spain: a Year-Long Longitudinal Study. *Appl Environ Microb* 79, 449-459, 2013.